

Product Data Sheet

Anti-CA242 Antibody (Detection)

Catalog #	Source	Reactivity	Applications
CAP1013-D	Mouse	H	ELISA, ICA, CLIA
Description	Mouse monoclonal detection antibody to CA242		
Immunogen	Recombinant protein of human CA242		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human CA242 protein. Does not cross-react with CA125, CEA, CA199 and CA724.		
Clonality	Monoclonal		
Conjugation			
Form	IgG1, liquid in PBS, pH 7.4.		
Dilution			
Gene Symbol			
Alternative Names			
Entrez Gene			
SwissProt			
Directions for Use	Recommended pair for sandwich immunoassay - Capture: CAP1013-C		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 2 years. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com