

Product Data Sheet

Anti-IGFBP1 Antibody (Detection)

Catalog #	Source	Reactivity	Applications
CAP1030-D	Mouse	H	ELISA, ICA, CLIA
Description	Mouse monoclonal detection antibody to IGFBP1		
Immunogen	Recombinant protein of human IGFBP1		
Purification	The antibody was purified by immunogen affinity chromatography.		
Specificity	Recognizes native and recombinant human IGFBP1 protein.		
Clonality	Monoclonal		
Conjugation			
Form	IgG1, liquid in PBS, pH 7.4, 0.05% Proclin 300.		
Dilution			
Gene Symbol	IGFBP1		
Alternative Names	IBP1; Insulin-like growth factor-binding protein 1; IBP-1; IGF-binding protein 1; IGFBP-1; Placental protein 12; PP12		
Entrez Gene	3484		
SwissProt	P08833		
Directions for Use	Recommended pair for sandwich immunoassay - Capture: CAP1030-C		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 2 years. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com