

Product Data Sheet

Anti-IL-5 Antibody (Capture)

Catalog #	Source	Reactivity	Applications
CAP1039-C	Rat	H	ELISA, ICA, CLIA
Description	Rat monoclonal capture antibody to IL-5		
Immunogen	Recombinant protein of human IL-5		
Purification	The antibody was purified by immunogen affinity chromatography.		
Specificity	Recognizes native and recombinant human IL-5 protein.		
Clonality	Monoclonal		
Conjugation			
Form	IgG1, liquid in PBS, pH 7.4, 0.05% Proclin 300.		
Dilution			
Gene Symbol	IL5		
Alternative Names	Interleukin-5; IL-5; B-cell differentiation factor I; Eosinophil differentiation factor; T-cell replacing factor; TRF		
Entrez Gene	3567		
SwissProt	P05113		
Directions for Use	Recommended pair for sandwich immunoassay - Detection: CAP1039-D		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 2 years. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com