

Product Data Sheet

Anti-sFlt-1 Antibody (Detection)

Catalog #	Source	Reactivity	Applications
CAP1047-D	Mouse	H	ELISA, ICA, CLIA
Description	Mouse monoclonal detection antibody to sFlt-1		
Immunogen	Recombinant protein of human sFlt-1		
Purification	The antibody was purified by immunogen affinity chromatography.		
Specificity	Recognizes native and recombinant human sFlt-1 protein.		
Clonality	Monoclonal		
Conjugation			
Form	IgG1, liquid in PBS, pH 7.4, 0.05% Proclin 300.		
Dilution			
Gene Symbol	FLT1		
Alternative Names	FLT; FRT; VEGFR1; Vascular endothelial growth factor receptor 1; VEGFR-1; Fms-like tyrosine kinase 1; FLT-1; Tyrosine-protein kinase FRT; Tyrosine-protein kinase receptor FLT; FLT; Vascular permeability factor receptor		
Entrez Gene	2321		
SwissProt	P17948		
Directions for Use	Recommended pair for sandwich immunoassay - Capture: CAP1047-C		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 2 years. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com