

Product Data Sheet

Anti-CXCL12 Antibody (Capture)

Catalog #	Source	Reactivity	Applications
CAP1071-C	Mouse	Human	SE
Description	Mouse monoclonal capture antibody to CXCL12		
Immunogen	Recombinant protein of human CXCL12		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human CXCL12 protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG2a. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	CXCL12		
Alternative Names	SDF1; SDF1A; SDF1B; Stromal cell-derived factor 1; SDF-1; hSDF-1; C-X-C motif chemokine 12; Intercrine reduced in hepatomas; IRH; hIRH; Pre-B cell growth-stimulating factor; PBSF		
Entrez Gene	6387		
SwissProt	P48061		
Directions for Use	Recommended pair for sandwich immunoassay - Detection: CAP1071-D		
Storage/Stability	Shipped at 4° C. Upon delivery aliquot and store at -20° C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com