

# Product Data Sheet

## AKT Blocking Peptide

| Catalog #                 | Source  | Reactivity | Applications |
|---------------------------|---|------------|--------------|
| CBP4560                   | Synthetic   | H, M, R, B | BL           |
| <b>Description</b>        | The peptide is used to block Anti-AKT Antibody (#CPA4560) reactivity.   |            |              |
| <b>Form</b>               | Lyophilized powder  |            |              |
| <b>Gene Symbol</b>        | AKT1; AKT2; AKT3  |            |              |
| <b>Alternative Names</b>  | PKB; RAC; RAC-alpha serine/threonine-protein kinase; Protein kinase B; PKB; Protein kinase B alpha; PKB alpha; Proto-oncogene c-Akt; RAC-PK-alpha                                 |            |              |
| <b>Entrez Gene</b>        | 207; 208; 10000 (Human); 11651; 11652; 23797 (Mouse); 24185; 25233; 29414 (Rat)   |            |              |
| <b>SwissProt</b>          | P31749; P31751; Q9Y243 (Human); P31750; Q60823; Q9WUA6 (Mouse); P47196; P47197; Q63484 (Rat)  |            |              |
| <b>Purity</b>             | >85%  |            |              |
| <b>Quality Control</b>    | The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.   |            |              |
| <b>Directions for Use</b> | Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |            |              |
| <b>Storage/Stability</b>  | Shipped at 4°C. Store at -20°C for one year.  |            |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)