

Product Data Sheet

**[Lys1015, 1024]-Thrombospondin-1 (1015-1024)
(Human, Bovine, Mouse)**

Catalog #	Source	Reactivity	Applications
CCP1791	Synthetic		
Description	Peptide to [Lys1015, 1024]-Thrombospondin-1 (1015-1024) (Human, Bovine, Mouse)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C68H105N17O12S		
Molecular Weight	1384.77		
Purity	> 95%		
Chemical Structure	Lys - Arg - Phe - Tyr - Val - Val - Met - Trp - Lys - Lys		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com