

Product Data Sheet

Chorionic Gonadotropin-Beta (109-145) (Human)

Catalog #	Source	Reactivity	Applications
CCP2208	Synthetic		
Description	Peptide to Chorionic Gonadotropin-Beta (109-145) (Human)		
Form	Lyophilized powder		
CAS Number	72979-70-5		
Molecular Formula	C167H264N46O58S1		
Molecular Weight	3876.22		
Purity	> 95%		
Chemical Structure	Thr - Cys - Asp - Asp - Pro - Arg - Phe - Gln - Asp - Ser - Ser - Ser - Ser - Lys - Ala - Pro - Pro - Pro - Ser - Leu - Pro - Ser - Pro - Ser - Arg - Leu - Pro - Gly - Pro - Ser - Asp - Thr - Pro - Ile - Leu - Pro - Gln		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com