

Product Data Sheet

IL-1b (208-240) (Human)

Catalog #	Source	Reactivity	Applications
CCP2505	Synthetic		
Description	Peptide to IL-1b (208-240) (Human)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C191H292N48O51S		
Molecular Weight	4108.81		
Purity	> 95%		
Chemical Structure	Lys - Lys - Lys - Met - Glu - Lys - Arg - Phe - Val - Phe - Asn - Lys - Ile - Glu - Ile - Asn - Asn - Lys - Leu - Glu - Phe - Glu - Ser - Ala - Gln - Phe - Pro - Asn - Trp - Tyr - Ile - Ser - Thr		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com