

PACAP-38 (28-38) (Human, Chicken, Mouse, Ovine, Porcine, Rat)

Catalog #	Source	Reactivity	Applications
CCP2715	Synthetic		

Description Peptide to PACAP-38 (28-38) (Human, Chicken, Mouse, Ovine, Porcine, Rat)

Form Lyophilized powder

CAS Number 160489-86-1

Molecular Formula C₆₁H₁₁₀N₂₄O₁₄


Molecular Weight 1403.71

Purity > 95%

Chemical Structure Gly - Lys - Arg - Tyr - Lys - Gln - Arg - Val - Lys - Asn - Lys - NH₂

Storage/Stability Shipped at 4°C. Store at -20°C for one year.

Structure



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com