

Product Data Sheet

Pancreatic Polypeptide (Rana temporaria)

Catalog #	Source	Reactivity	Applications
CCP2724	Synthetic		
Description	Peptide to Pancreatic Polypeptide (Rana temporaria)		
Form	Lyophilized powder		
CAS Number	132187-74-7		
Molecular Formula	C192H276N52O58		
Molecular Weight	4240.65		
Purity	> 95%		
Chemical Structure	Ala - Pro - Ser - Glu - Pro - His - His - Pro - Gly - Asp - Gln - Ala - Thr - Gln - Asp - Gln - Leu - Ala - Gln - Tyr - Tyr - Ser - Asp - Leu - Tyr - Gln - Tyr - Ile - Thr - Phe - Val - Thr - Arg - Pro - Arg - Phe - NH ₂		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com