

Product Data Sheet

Pannexin-1 Fragment (4512)

Catalog #	Source	Reactivity	Applications
CCP2725	Synthetic		
Description		Peptide to Pannexin-1 Fragment (4512)	
Form		Lyophilized powder	
CAS Number			
Molecular Formula		C55H94N18O24	
Molecular Weight		1391.47	
Purity		> 95%	
Chemical Structure		Val - Gln - Gln - Lys - Ser - Ser - Leu - Gln - Ser - Glu - Ser - Gly - Asn	
Storage/Stability		Shipped at 4°C. Store at -20°C for one year.	
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com