


Peptide YY (3-36) (Canine, Mouse, Porcine, Rat)

Catalog #	Source	Reactivity	Applications
CCP2761	Synthetic		
Description	Peptide to Peptide YY (3-36) (Canine, Mouse, Porcine, Rat)		
Form	Lyophilized powder		
CAS Number	126339-09-1		
Molecular Formula	C176H272N52O54		
Molecular Weight	3980.35		
Purity	> 95%		
Chemical Structure	Ala - Lys - Pro - Glu - Ala - Pro - Gly - Glu - Asp - Ala - Ser - Pro - Glu - Glu - Leu - Ser - Arg - Tyr - Tyr - Ala - Ser - Leu - Arg - His - Tyr - Leu - Asn - Leu - Val - Thr - Arg - Gln - Arg - Tyr - NH ₂		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Structure



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com