

Product Data Sheet

TIP-39

Catalog #	Source	Reactivity	Applications
CCP2937	Synthetic		
Description	Peptide to TIP-39		
Form	Lyophilized powder		
CAS Number	277302-47-3		
Molecular Formula	C202H325N61O54S		
Molecular Weight	4504.28		
Purity	> 95%		
Chemical Structure	Ser - Leu - Ala - Leu - Ala - Asp - Asp - Ala - Ala - Phe - Arg - Glu - Arg - Ala - Arg - Leu - Leu - Ala - Ala - Leu - Glu - Arg - Arg - His - Trp - Leu - Asn - Ser - Tyr - Met - His - Lys - Leu - Leu - Val - Leu - Asp - Ala - Pro		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com