

VIP-Lys (Biotin) (Human, Porcine, Rat)

| Catalog # | Source | Reactivity | Applications |
|---------------------------|---|------------|--------------|
| CCP2985 | Synthetic | | |
| Description | Peptide to VIP-Lys (Biotin) (Human, Porcine, Rat) | | |
| Form | Lyophilized powder | | |
| CAS Number | | | |
| Molecular Formula | C163H263N47O46S2 | | |
| Molecular Weight | 3681.33 | | |
| Purity | > 95% | | |
| Chemical Structure | His - Ser - Asp - Ala - Val - Phe - Thr - Asp - Asn - Tyr - Thr - Arg - Leu - Arg - Lys - Gln - Met - Ala - Val - Lys - Lys - Tyr - Leu - Asn - Ser - Ile - Leu - Asn - Lys (Biotin) | | |
| Storage/Stability | Shipped at 4°C. Store at -20°C for one year. | | |
| Structure | | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com