

Product Data Sheet

Tyr-CRF (Human, Rat)

Catalog #	Source	Reactivity	Applications
CCP3002	Synthetic		
Description	Peptide to Tyr-CRF (Human, Rat)		
Form	Lyophilized powder		
CAS Number	100513-58-4		
Molecular Formula	C217H353N61O65S2		
Molecular Weight	4920.73		
Purity	> 95%		
Chemical Structure	Tyr - Ser - Glu - Glu - Pro - Pro - Ile - Ser - Leu - Asp - Leu - Thr - Phe - His - Leu - Leu - Arg - Glu - Val - Leu - Glu - Met - Ala - Arg - Ala - Glu - Gln - Leu - Ala - Gln - Gln - Ala - His - Ser - Asn - Arg - Lys - Leu - Met - Glu - Ile - Ile - NH2		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com