

Product Data Sheet

[Tyr0]-Stresscopin-Related Peptide (Human)

Catalog #	Source	Reactivity	Applications
CCP3008	Synthetic		
Description	Peptide to [Tyr0]-Stresscopin-Related Peptide (Human)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C214H367N69O59		
Molecular Weight	4850.74		
Purity	> 95%		
Chemical Structure	Tyr - His - Pro - Gly - Ser - Arg - Ile - Val - Leu - Ser - Leu - Asp - Val - Pro - Ile - Gly - Leu - Leu - Gln - Ile - Leu - Leu - Glu - Gln - Ala - Arg - Ala - Arg - Ala - Ala - Arg - Glu - Gln - Ala - Thr - Thr - Asn - Ala - Arg - Ile - Leu - Ala - Arg		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com