

Product Data Sheet

Nesfatin-1 (Mouse)

Catalog #	Source	Reactivity	Applications
CCP3045	Synthetic		
Description	Peptide to Nesfatin-1 (Mouse)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C424H683N117O137		
Molecular Weight	9611.87		
Purity	> 95%		
Chemical Structure	Val - Pro - Ile - Asp - Val - Asp - Lys - Thr - Lys - Val - His - Asn - Thr - Glu - Pro - Val - Glu - Asn - Ala - Arg - Ile - Glu - Pro - Pro - Asp - Thr - Gly - Leu - Tyr - Tyr - Asp - Glu - Tyr - Leu - Lys - Gln - Val - Ile - Glu - Val - Leu - Glu - Thr		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com