


Sarafotoxin S6b

Catalog #	Source	Reactivity	Applications
CCP3049	Synthetic		
Description	Peptide to Sarafotoxin S6b		
Form	Lyophilized powder		
CAS Number	120972-53-4		
Molecular Formula	C ₁₁₀ H ₁₅₉ N ₂₇ O ₃₄ S ₅		
Molecular Weight	2563.96		
Purity	> 95%		
Chemical Structure	Cys - Ser - Cys - Lys - Asp - Met - Thr - Asp - Lys - Glu - Cys - Leu - Tyr - Phe - Cys - His - Gln - Asp - Val - Ile - Trp (Disulfide bridge Cys1 - Cys15, Cys3 - Cys11)		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com