

Product Data Sheet

hsa-miR-450a-2-3p miRNA Mimic

Catalog #	Source	Reactivity	Applications
CMH2564	Synthetic	H	
Description	Synthetic miRNA mimics are used to overexpress hsa-miR-450a-2-3p by transfection.		
Specificity	<p>miRNA Mimics are chemically modified, double-stranded miRNA-like RNA which are designed to copy the functionality of mature endogenous miRNA upon transfection. At the 5'-end, it is synthesized with a partially complementary motif to 3'UTR end of the target gene, which allows the miRNA mimic to specifically bind to the target. Transfection of mimics followed by downstream gene expression analysis or phenotypic analysis, is performed to elucidate the targets and roles of particular miRNAs.</p>		
Form	Lyophilized powder		
Gene Symbol	hsa-miR-450a-2-3p		
Accession No.	MIMAT0031074		
Components	This synthetic miRNA is based on the mature miRNA sequence. It does not contain the full precursor miRNA stem-loop.		
Directions for Use	We recommend re-suspending the lyophilized synthetic miRNA using DNase and RNase-free ddH ₂ O. To make a 100 uM stock, dissolve the lyophilized powder using 50 ul of ddH ₂ O.		
Storage/Stability	Shipped at 4 °C. Store at -20 °C for one year. Avoid freeze-thaw cycles after reconstitution.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com