

Product Data Sheet

mmu-miR-1970 miRNA Mimic

Catalog #	Source	Reactivity	Applications
CMM0652	Synthetic	M	
Description	Synthetic miRNA mimics are used to overexpress mmu-miR-1970 by transfection.		
Specificity	<p>miRNA Mimics are chemically modified, double-stranded miRNA-like RNA which are designed to copy the functionality of mature endogenous miRNA upon transfection. At the 5'-end, it is synthesized with a partially complementary motif to 3'UTR end of the target gene, which allows the miRNA mimic to specifically bind to the target. Transfection of mimics followed by downstream gene expression analysis or phenotypic analysis, is performed to elucidate the targets and roles of particular miRNAs.</p>		
Form	Lyophilized powder		
Gene Symbol	mmu-miR-1970		
Accession No.	MIMAT0009444		
Components	This synthetic miRNA is based on the mature miRNA sequence. It does not contain the full precursor miRNA stem-loop.		
Directions for Use	We recommend re-suspending the lyophilized synthetic miRNA using DNase and RNase-free ddH ₂ O. To make a 100 uM stock, dissolve the lyophilized powder using 50 ul of ddH ₂ O.		
Storage/Stability	Shipped at 4 °C. Store at -20 °C for one year. Avoid freeze-thaw cycles after reconstitution.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com