

## L-methionine gamma-lyase

Catalog #	Source	Reactivity	Applications
CZM1116	E.coli		
<b>Description</b>	Recombinant L-methionine gamma-lyase from E.coli		
<b>Specificity</b>	>40 U/mg protein, >15 U/mg powder		
<b>Form</b>	Light yellow lyophilized powder		
<b>Alternative Names</b>	EC 4.4.1.11; L-methioninase; methionine lyase; methioninase; methionine dethiomethylase; L-methionine $\gamma$ -lyase; L-methionine methanethiol-lyase (deaminating)		
<b>CAS Number</b>	42616-25-1		
<b>Molecular Weight</b>			
<b>Purity</b>	>90%		
<b>Unit Definition</b>	One unit is defined as the amount of enzyme will release 1 micromole alpha-ketobutyrate from L-methionine per minute at pH 8.0 and 37°C.		
<b>Storage/Stability</b>	Shipped at 4°C. Store at -20°C for 2 years.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)