

# Product Data Sheet

## Glucokinase

Catalog #	Source	Reactivity	Applications
CZM1122	E.coli		
<b>Description</b>	Recombinant Glucokinase from E.coli		
<b>Specificity</b>	>300 U/mg protein, >20 U/mg powder		
<b>Form</b>	White lyophilized powder		
<b>Alternative Names</b>	EC 2.7.1.2; GCK; FGQTL3; GK; GLK; HHF3; HK4; HKIV; HXKP; LGLK; MODY2		
<b>CAS Number</b>	9001-36-9		
<b>Molecular Weight</b>			
<b>Purity</b>	>90%		
<b>Unit Definition</b>	One unit is defined as the amount of enzyme will convert 1 nanomole NADP to NADPH per minute at 30°C.		
<b>Storage/Stability</b>	Shipped at 4°C. Store at -20°C for 2 years.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com